

Wykeham Solo

Instruction Manual

Polygiene

All the sections of the mesh on the Wykeham Solo are coated in Polygiene to ensure there are no health and hygiene risks. Polygiene lasts the lifetime of the fabric by being permanent odour control, inhibiting the growth of odour-causing bacteria.

The patented technology helps to protect the mesh against microbes by building one of nature's own antimicrobial agents into various types of materials & products. Harmful microbes deposited onto these surfaces become incapable of replicating themselves through a process known as:

"The Polygiene Effect"

Polygiene technology is based on the natural ability of silver salt to inhibit the growth of a wide range of harmful microbes. In products made with Polygiene, silver ions continuously migrate to the surface where the silver ions effectively suppress microbial growth. This activity provides protection 24 hours a day, 7 days a week throughout the normal lifespan of the product.

www.easycare systems.co.uk

easyCare
systems

Special products and Service



The Wykeham Solo is designed to aid bathing and showering. In conjunction with any of the Daily Care cradles the base allows easy room to room transportation.

Before operation ensure:

- All 4 castor wheels are firmly attached to the base with no damage
- The mesh is clean, undamaged and correctly attached to the frame via Velcro and zips.

Operation

When hoisting a person onto or off the Wykeham unit ensure all 4 brakes are applied to the castor wheels before hand. If the Wykeham unit is stationary Daily Care advise to apply the brakes to the castor wheels.

The Wykeham unit is made from Stainless Steel and is powder coated, this allows the product to be placed directly into a shower, wet room and the cradle into a bath. With the design of the Wykeham cradle being relatively flat, the cradle can be placed directly into a bath without the need for any accessories as the two support bars act at the stand.

When placing the cradle into a bath firstly check the castor wheel brakes are applied. The Red Triangles can then be rotated so the cradle can be removed from the base. Lift up the cradle and carefully place it into a bath. Once the cradle is correctly and securely placed in the bath, with the use of a hoist and sling the user can be carefully lowered onto the cradle.

Note: do not attempt to hoist or attach a hoist or hoist straps directly to the Wykeham unit / cradle.

Mesh

The Wykeham Cradle can be detached from the Base via the two rotating Red Triangles, these rotate allowing the cradle to be lifted off the Base.

It is advised that the cradle mesh is not changed when there is someone on the cradle. The simplest way of changing the mesh is to detach the cradle from the Base, spin it upside down and place it against a sturdy surface. This gives access to underneath the cradle where the mesh Velcro and zips can be accessed. These Velcro and zips can be undone for washing of the mesh, when washing please follow the instructions on the mesh.





Special Products and Service

Contact Us

Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk