



## MANGAR Eagle Lifting Cushion

The Eagle lifting cushion will both sit up and lift a fallen person. It has a small footprint and can be used in small spaces such as a bathroom.

### EASY TO FIT

- Lifting cushion is inflated using an Airflo Compressor with a simple hand control
- Simple to use, requiring very little training
- Suitable for indoor and outdoor use
- Can be used independently or assisted

### SIMPLE TO MAINTAIN

- Clean using a proprietary non-abrasive mild liquid cleaner
- Easy to deflate, fold away and store in the carry bag provided
- Does not need a LOLER certificate
- The material is designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

### SAFE TO USE

- Lifts the frailest person up to 35 stone/ 220 kg from the floor
- Exceptionally rigid and supportive to back
- Ideal for use in supported living apartments and home care

### IMPROVES WELL-BEING

- Provides a safe, dignified lift for a fallen person
- Reduces the risk of musculoskeletal injury for the carer
- The backrest gives additional support to raise the person from lying to sitting
- Small footprint is ideal for use in a home environment

[www.easycare systems.co.uk](http://www.easycare systems.co.uk)

### TECHNICAL SPECIFICATIONS

Max user weight	220 kg
-----------------	--------

Max seat height	56 cm
-----------------	-------

Min height	Flat
------------	------

Seat depth	57 cm
------------	-------

Total length	112 cm
--------------	--------

Width	64 cm
-------	-------

Total product weight	4.5 kg
----------------------	--------

Back/headrest length	52 cm
----------------------	-------

Airflo Compressor	24
-------------------	----

Lifting performance per charge	up to 6 lifts
--------------------------------	---------------



## Features

- 1. Easy to setup** - connect the four hoses from the hand control to the corresponding numbers and coloured connectors on the Eagle.
- 2. The backrest** provides support to the back and head.
- 3. Made from very strong polyurethane coated nylon material.**
- 4. Each layer inflates individually** ideal for single-handed care
- 5. Sides of the Eagle can be folded up** to allow manoeuvring into tight spaces.



## Accessories

- Slidesheet
- Eagle Carry Bag



## For more information visit

Units 5 and 6 Lower Vicarage  
Road Southampton Hampshire  
SO19 7RJ  
T: 02380 420678  
E: info@easycare systems.co.uk

easyCare  
systems

[www.easycare systems.co.uk](http://www.easycare systems.co.uk)