

Moray Combo™


Our ranges of Dynamic Hybrids provide the latest in innovation and technology, both products are dual mode comprising of Alternating and Constant Low Pressure, a range of comfort settings eliminate the need to set to Patients weight, a transport mode allows the mattress to remain inflated when not connected to the Control unit and a Static mode which ensures a firm surface for Patient care returning to the previous cycle after 30 mins, our Visco heel-zone has been created to ensure minimum shear and friction takes place on the most vulnerable parts of the body.

Our Advanced version provides an inner surface of High density Visco which covers the Dynamic Cells ensuring maximum patient comfort whilst receiving the benefits of both alternating and CLP therapeutic surfaces.

**If you have any questions we are here to help.
Call us on 0800 193 6030.**

PRODUCT FEATURES

- Alternating and Constant Low Pressure
- Transport mode allows the mattress to remain inflated when not connected to the Control unit
- Visco heel-zone ensures minimum shear and friction on the most vulnerable parts of the body



Features and Benefits

- Fully Welded, Multi stretch- vapour permeable fabric allows high moisture vapour transfer rate (MVTR), maximising airflow to patient skin surface.
- Leading Design/Safer Heel Solution. The clinically proven lower leg section, redistributes pressure from the heel to the upper leg (calve).
- CME Foam surround and base maintains pelvic stability for transfer and aids moving and handling processes.
- High Quality Foam Filled Cell Technology. Foam filled, alternating cells that automatically inflate to a comfort level in line with the very best pressure relieving surfaces.
- Active Therapy Mode (WITH PUMP). A digital pump can also be fitted to the mattress, which in turn provides dynamic or alternating pressure therapy.
- Comfort Mode Selection. The pump is designed to offer ultra-low-pressure delivery, for light weight, frail individuals, including Paediatric patients.

TECHNICAL INFORMATION

Product Code	HER-COMBO-777 For dual mode with CLP use code HER-COMBO-777-CLP
Risk Level	Very High Risk
Maximum Patient Weight	35 stone 222kg
Mattress Size	198 x 88 x 16cm
Active Therapy mode (with pump)	
Comfort mode selection	
High quality foam filled cell technology	
CME foam surround & base	
Unique CPR function/Action	

WIN  CARE

easy  care
systems

Special Products and Service

Contact Us
Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk