

MESH SIDE RAILS FOR ULTRA-LOW BED



Product Code: HLDWUL110

Description: Mesh Side Rails for Ultra-Low Bed

- Intended to be used to provide security to the occupant of the bed while the bed is raised to height for giving care.
- The Woburn Ultra-Low bed is designed to be lowered to the floor when the occupant in the bed is left unattended, this is to reduce the risk of a fall from the bed or potential entrapment from a side rail system.
- When the bed is lowered to the floor it is recommended the side rail system is locked down by connecting the two halves of the plastic buckle at the centre on the outside edge.
- It can be used with the enclosed head and foot boards with some adjustments.

Note:

- This system is only compatible with the Harvest Woburn Ultra Low bed
- Carry out a risk assessment to ensure the bumper is suitable for the users' requirements
- In the event of the bumper system becoming damaged discontinue use immediately and contact Harvest Healthcare for advice.