

How to fit your stick clip (HX15153R_1)

Keep your reacher/grabber with you wherever you go.



This stick clip is specifically designed to fit snugly on all 'Handi-Grip' & 'ArthriGrip' reachers

To fit on to your reacher:



rabbit ears

**Pinch the 'rabbit ears'
to open the clip**



**Push onto the side
of the reacher bar**



**Make sure it 'clunk-
clicks' into place**

Clip it on to a crutch, chair leg
walking pole/stick/frame or
your wheelchair – simple.

For more information
and other useful products
please visit:
www.easycare systems.co.uk

