


# How to fit your stick clip (HX15151R\_1)

Keep your reacher/grabber with you wherever you go.

**This stick clip is specifically designed to fit snugly on all 'Classic' reachers**

**To fit on to your reacher:**



**Pinch the 'rabbit ears'  
to open the clip**



**Push onto the side  
of the reacher bar**



**Make sure it 'clunk-  
clicks' into place**

*Clip it on to a crutch, chair leg walking pole/stick/frame or your wheelchair – simple.*

*For more information  
and other useful products  
please visit:  
[www.easycare systems.co.uk](http://www.easycare systems.co.uk)*

Easy Care Systems  
Units 5 and 6 Lower Vicarage Road  
Southampton  
Hampshire  
SO19 7RJ  
T: 02380 420678  
E: [info@easycare systems.co.uk](mailto:info@easycare systems.co.uk)

