


Handi-Grip[®] MAX



Operating instructions

Hold jaw and collar, and pull



Rotate jaw 90° then push back

HA8326 26"/66cm Handi-Grip Max HA8332 32"/81cm Handi-Grip Max

Gripcert gel handle

Lock button

To close; pull back lock and trigger together

To release; pull trigger only

CAUTION

- GB - Don't use for picking up items above 60°C
- Do not use as walking aid
- Do not lift objects greater than 1kg




HX80083GB_2

Handi-Grip[®] MAX



Operating instructions

Hold jaw and collar, and pull



Rotate jaw 90° then push back

HA8326 26"/66cm Handi-Grip Max HA8332 32"/81cm Handi-Grip Max

Gripcert gel handle

Lock button

To close; pull back lock and trigger together

To release; pull trigger only

CAUTION

- GB - Don't use for picking up items above 60°C
- Do not use as walking aid
- Do not lift objects greater than 1kg



HX80083GB_2

Easy Care Systems
Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk

www.easycare systems.co.uk



easyCare systems

Helping Hand Reablement
A division of The Helping Hand Company

Easy Care Systems
Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk

www.easycare systems.co.uk



easyCare systems

Helping Hand Reablement
A division of The Helping Hand Company