easystems Shower Chairs, Stretchers and Trolleys



www.easycaresystems.co.uk

QUALITY. D

Backrests and handles

All models are available with either horizontal push handles or vertical handles for use in restricted areas or over WCs. The backrest is pliable and moulds to the user while providing excellent support.

FREECOTE

All Freeway chairs are treated with Freecote, our unique antibacterial and anti-microbial coating. The coating is durable and prolongs the life of the chair.

Frames

Available in three sizes:

- 'A' 540mm (21") wide
- 'B' 490mm (19") wide
- 'C' 440mm (17") wide

High performance wheels

With stainless steel bearings and puncture proof tyres.

Armrests

With soft but supportive tops. The armrests detach or swing away for easy access and transfer.

Seats

Fully padded, horseshoe or full ring seats with different sized apertures that are interchangeable.

Footrests

They flip up and are height adjustable. They can also be easily removed or swung conveniently away from the user's foot area.

Brakes

Self-propelled models have large wheel brakes fitted as standard.

Castors Highest quality large diameter plastic castors with stainless steel bearings.

URABILITY. HYGIENE.

Bespoke, cost effective toileting and showering solutions

Manufactured in the UK from only the highest quality materials, the comprehensive range of Freeway assisted or self-propelled shower, toilet and commode chairs are designed to provide the most effective solution to a client's toileting and showering needs.

Freeway Shower Chairs are of **modular design**, and can be configured to a choice of interchangeable seats, backrests, armrests and footrests. A multitude of accessories are available to ensure the most bespoke solution.

All adult models within the range are specifically designed to fit over a standard toilet or commode, with three widths available; 'A' - 540mm (21"), 'B' - 490mm (19") and 'C' - 440mm (17"). Paediatric, Bariatric and models to fit over automatic toilets are also included in the range.

The interchangeable makeup of Freeway Shower Chairs makes the range one of the most cost effective on the market, as accessories and parts can be added to suit the changing needs of a client. Frames can be lowered/ raised; rear wheels set back/forward and the frame depth extended/shortened.

Where there are special needs, Freeway can manufacture a bespoke shower chair subject to an individual assessment.

All standard models are weight tested to 200kg (30 stone)* and come with a standard 12 month warranty.

* excluding the Freeway T80 which is tested to 160kg (25 stone)







FREEWAY Shower



Freeway T20 **Moulded Shower** Seat Chair

Available as a static model or with four castors. It has a moulded polypropylene seat and fixed arms.



Freeway T30 **Free Standing Shower Chair**

The Freeway T30 is a static shower chair with slip resistant rubber feet or floor fixing options available.



Freeway T40 Assistant **Propelled Shower Chair**

Versatility is the key with the T40 as it can be raised/lowered, extended forward, rear castors set back for amputee use, to name but a few of its applications. Features four total stop castors and footrests as standard.



The Freeway T40 **COM** is a version of the standard T40 which has been modified to work seamlessly with the Clos-o-Mat Samoa toilet.

The Freeway T40 AUTO is for use with automatic toilets including the Geberit 8000 range and Clos-o-Mat Palma Vita. Features pram handle backrest and seat with skirt as standard.



Chairs

All models available as shower, toilet or commode chairs. (Excluding T20)



Freeway T50 Front Wheel Self Propelled Shower Chair

The Freeway T50 features 510mm (20") braked front wheels and includes footrests as standard. It is ideal for operating in tight spaces due to the positioning of the castors.

Freeway T60 Self Propelled Shower Chair

Benefits from a short overall length making it ideal for small shower cubicles. The chair's low braked wheel facilitates seamless side transfers as the top of the wheels are below the seat height of the chair. Anti-tip bars are featured as standard. Suitable for self-propelled or assistant propelled as the 510mm (20") rear wheel aids manoeuvring over a lip or small step.





Freeway T70 Self Propelled Shower Chair

Benefits from large braked rear wheels making this chair easy to self-propel.

The standard self-propelled model with 560mm (22") rear wheels has easily removed push rims and footrests fitted.

Freeway T70 COM: Modified to work seamlessly with the Clos-o-Mat Samoa toilet.

Freeway T70 AUTO: Designed specifically for use with automatic toilets including the Geberit 8000 range and Clos-o-Mat Palma Vita.

Freeway T70 Side Access: Equipped with left/right side opening and adapted, swing away armrests.

Reclining







Freeway T80 Reclining Shower Chair

The Freeway T80 is highly versatile and the reclined position has significant advantages for a wide range of conditions, increasing comfort and relieving pressure.

Such conditions include:

- Multiple sclerosis
- Spinal injuries
- Motor Neurone Disease
- Stroke
- Celebral palsy
- Other neurological conditions

The chair reclines through an infinite number of positions up to 45 degrees. The tilt action is easy and light to operate.

The Freeway T80 is available as assistant propelled or self propelled and has a Safe Working Load of 160kg (25 stone).

An extensive range of accessories are available for this model such as a full or half height detachable headrest.

Freeway T80 COM

A variation of the T80 reclining shower chair, designed to work seamlessly with the Clos-o-Mat Samoa toilet.

Freeway T80 AUTO

A version of the standard T80 which has been modified for use with automatic toilets including the Geberit 8000 range and Closo-Mat Palma Vita.

Paediatric

The Freeway Paediatric Shower Chair range features a seat width of 360mm and a seat depth of 325mm

Freeway T80 Paediatric Reclining Shower Chair

The Freeway T80 Paediatric is a reclining shower or toilet/commode chair that has been purpose built for the needs of younger clients.

Freeway T90 Paediatric Shower Chair

The Freeway T90 Shower Chair/Commode features a seat height of 460mm. Footplates are flip up, removable and height/depth adjustable. The T90 is available as an assisted or self propelled model.







The **Freeway T90 WC** is a variation of the above, but for usage as a toilet commode - with a seat height of 560mm, allowing clearance over standard toilets and automatic bidets. NB: Self propelled model is wider due to the splayed frame

All models available as shower, toilet or commode chairs.

Bariatric

Freeway T100 Bariatric Shower Chair



The Freeway T100 carries up to 350kgs (55 stone) and is suitable for a variety of applications for bariatric clients. Seat dimensions are 790mm wide and 550mm deep.

The T100 benefits from increased diameter tubing, special heavy duty castors, and heavy duty footrests with heel supports which are height adjustable, removable and swing away. The T100 is also available in smaller width frames than the standard.

Adaptations and A

Freeway Shower Chairs can have special adaptations, accessories and alterations to suit specific needs. Some of the options shown here may be fitted to an existing chair, others will be available only on a new build. Many will need adapting to your special size/situation requirements, subject to an individual assessment.

Armrest Options





Raised/lowered armrests

Flip up arms (only on new chair orders)

Backrest Options



Gull wing arms (only on new chair orders)



Extended forward armrests



Gutter armrests (with or without strap supports)



Padded style armrests



5° reclined backrest

Footrest Options



Extended backrest

Detachable headrest




Clip on padding for backrests and armrests



Curved padded head support



Raised footrests (60, 90



Dartex footplate sleeve



Heavy duty footrests



ABS foot tray



support



Integral heel supports with Clip on padded foot tray Velcro toe/foot straps



Padded ABS foot tray





Individual clip on foot tray



Adjustable leg support and footrest




Wrap around vinyl footrest covers



Adjustable angle footrest



Adjustable stump support



8

ccessories 7

Straps, Harnesses and Positional Aids



Waist band Velcro fastening



Waist band buckle fastening



Full harness



Ankle strap



Vinyl calf support



Armrest Padding



Adjustable lateral Supports



Thoracic wedge

Seat Options

٠

٠

٠

٠

٠



Trunk supports, padded



Trunk supports unpadded

Hygiene



Knee block (only on new chairs)



Black curved thoracic support (only with extended backrests)



Curved thoracic supports (can be purchased as a retro fit kit for all T80 length backrests)





Cup shaped splashguard



Lower commode pan tray



Pommel (photo shows adapted for use on horseshoe seat)

The following are examples of this type of seat.

Padding: Increased or extra comfort padding

• Dartex upholstery (breathable four way stretch fabric allows more significant moulding of seat to user

any Freeway Shower Chair.

Reduced padding on seats

and may help with pressure relief) Different sized/shaped apertures

Wedged padding

Commode pan

A padded seat is a major contributor to the comfort and support offered by

The huge choice of standard seats means that there is a seat for most situations yet we frequently supply custom seats for special situations.





Other Options



Counter balance weights for double amputee



Capstans on wheels

Extended brake levers



T30 floor fixing kit



Dartex Covered Seat




Extended seat depths (over 75mm involves frame alteration)








Build your own FREE



EWAY Shower Chair



STEP 6

TRANSPORT

Many Freeway Shower Chairs are available as a transportable option so that the frame can be easily dismantled without the need for tools, for ease of transport or storage.

For example:



T40 Assembled



T40 Disassembled (frame only)

Please speak to your representative for a full list of shower chairs with transportable option.

Transportable option





ACCESSORIES

Select your accessories

Please quote the product name of the accessories that you require:

Your Freeway Shower Chair is complete!

Pass your completed form to your local representative or call your local service centre.

Shower Chair Tecl

Below are the technical specifications of the most popular Freeway Shower Chair models. Please consult our website, www.prismmedical.co.uk for the full range of dimensions.

- 1 = Overall maximum height.
- 2 = Overall maximum width.
- 3 = Length including footrests (if fitted) at lowest position.
- 4 = Length including footrests (if fitted) at highest position.
- 5 = Width of standard seat.
- 6 = Height to top of seat.

- 7 = Height clearence.
- 8 = Seat depth (standard seat).
- 9 = Top of armrest to top of seat.
- 10 = Width clearance.
- 11 = Top of backrest to seat.
- 12 = Distance between armrests at rear.
- 13 = Distance between armrests at front.

тзо	12		95	9.		ि	P	Ð.	<u>.</u> 9?				
	1	2	3	4	5	6	7	8	9	10	11	12	13
T30A (SC only)	920	640	570	570	540	460	345	400	180	485	450	470	530
T30B (SC only)	920	590	570	570	490	460	345	400	180	435	450	420	480
T30C (SC only)	890	540	570	570	440	460	345	400	180	385	450	370	430
T40	1	2	3	4	5	<u>6</u>	ی ح ۲	8	* * 05	10		12	13
T40A	980	640	1050	960	540	560	445	400	180	485	450	470	530
T40B	980	590	1050	960	490	560	445	400	180	435	450	420	480
T40C	980	540	1050	960	440	560	445	400	180	385	450	370	430
T40A AUTO	1000	590	1050	960	540	560	480	400	180	485	450	460	530
T40B AUTO	1000	540	1050	960	490	560	480	400	180	435	450	410	480
T40C AUTO	1000	540	1050	960	440	560	480	400	180	385	450	410	430
T40A COM	980	690	1050	960	540	560	470	400	180	485	450	470	530
T40B COM	980	590	1050	960	490	560	470	400	180	435	450	420	480
T40C COM	980	540	1050	960	440	560	470	400	180	385	450	370	430
т50	192	2	3	4	5	6 6	<u>Б</u> . 7	8	*	10		12	13
T50A	980	730	940	900	540	560	445	400	180	485	450	470	530
T50B	980	680	940	900	490	560	445	400	180	435	450	420	480
T50C	980	630	940	900	440	560	445	400	180	385	450	370	430
	,		`	^							··		
т60	9	H	9.	9.	Ē	ਉਣਾ	P	£.	<u>.</u> 9?	Ģ		Ē	
-Sellen	1	2	3	4	5	6	7	8	9	10	11	12	13
T60A	980	710	1050	970	540	550	435	400	180	485	450	470	530
T60B	980	660	1050	970	490	550	435	400	180	435	450	420	480
T60C	980	610	1050	970	440	550	435	400	180	385	450	370	430

nnical

т70	B	<i>5</i> 7	Ð	9,	95		ਨਾ	P D	> B	» * }	2	Q			
	1		2	3	4	5	6	7	8	9		10	11	12	13
T70A	98	80	750	1150	1100	540	560	445	400	18	30 ·	485	450	470	530
T70B	98	0	700	1150	1100	490	560	445	400	18	30 ·	435	450	420	480
T70C	98	30	650	1150	1100	440	560	445	400	18	30 3	385	450	370	430
T70A AUTO	10	00	790	1140	1060	540	560	480	400	18	30 ·	485	450	460	530
T70B AUTO	10	00	740	1140	1060	490	560	480	400	18	30 ·	435	450	410	480
T70C AUTO	10	00	690	1140	1060	440	560	480	400	18	30 3	385	450	410	430
T70A COM	98	80	775	1110	1020	540	560	470	400	18	30 ·	485	450	470	530
T70B COM	98	80	725	1110	1020	490	560	470	400	18	30 ·	435	450	420	480
T70C COM	98	80	675	1110	1020	440	560	470	400	18	30 .	385	450	370	430
T90 Paediatric	13	<i>3</i> 7	H	<u>P</u>	95		रु	P.	× 8	× *	22				
	1		2	3	4	5	6	7	8	9		10	11	12	13
T90 Assistant Prop	92	20	620	730	N/A	360	460	N/A	330	13	1 08	N/A	440	300	340
T90 Self Propelled	92	20	620	880	N/A	360	460	N/A	330	13	1 08	N/A	440	300	340
T90WC Assistant Pr	. 98	30	470	680	N/A	360	560	470	330	13	50 ·	430	440	320	320
T90WC Self Prop.	98	80	680	880	N/A	360	560	470	330	13	50 ·	430	440	320	320
T100 Bariatric	13	<i>5</i> 7	H	<u>þ</u> ,	9.		97	j.	~ B	*	22	$\overline{\Omega}$			
Darraciae	1		2	3	4	5	6	7	8	9)	10	11	12	13
T100	99	0	900	1220	1170	790	570	450	550	17	70	730	430	720	760
T100FR	99	0	900	1220	1170	790	570	450	550	17	70 É	730	430	720	760
		i i		1								1			
T80 Baalinar	A	A	員	A	A	A	Ħ	科 1	A	A	Ą	ħ	, "員	員	員
Recliner	1	2	3	4 A	4 B	4C	5	6	7	8	9	10	11	12	13
T80A	1125	1440	640	1040	1300	1500	540	570	445	400	190	480	535	470	530
T80B	1125	1440	590	1040	1300	1500	490	570	445	400	190	430	535	420	480
T80C	1125	1440	540	1040	1300	1500	440	570	445	400	190	380	535	370	430
T80A COM	1125	1440		1040	1300	1500	540	570	470	400	190	480		460	530
T80B COM	1125	1440		1040	1300	1500	490	570	470	400	190	430		410	480
T80C COM	1125	1440		1040	1300	1500	440	570	470	400	190	380		410	430
T80A AUTO	1125	1440		1040	1300	1500	540	570	470	400	190	480	_	470	530
T80B AUTO	1125	1440		1040	1300	1500	490	570	470	400	190	430	_	420	480
T80C AUTO T80 Paediatric	1125 1050	1440 1345	+	1040 760	1300 1050	1500 1170	440 360	570 570	470 470	400	190 145	430 430		370 360	430
Touraculatric	1030	1343	390	700	1050	11/0	500	570	470	325	143	430	550	500	360

T80 Specification Key

1 = Height.

2 = Height with detachable headrest.

3 = Width.

4A = Length in upright position.

4B = Length with max. recline including footrests

at highest position.

4C = Length max. recline and detachable headrest fitted.

5 = Seat width.

6 = Height to top of seat.

7 = Height clearance.

8 = Seat depth (standard seat).

9 = Top of armrest to top of seat.

10 = Width clearance.

11 = Top of backrest to top of seat (without detachable

headrest fitted - with headrest add 300mm).

12 = Distance between armrests at rear.

13 = Distance between armrests at front.

DURABLE. MANOE

FREECOTE

All Freeway trolleys and stretchers are treated with Freecote, our unique anti-bacterial and antimicrobial coating. **Foot Pump** Height adjustable trolleys operated by easy to use foot pump.

> **Brakes** All trolleys have braked castors fitted as standard.

Castors Highest quality large diameter plastic castors with stainless steel bearings.

www.easycaresystems.co.uk

VREABLE. PRACTICAL.

Shower trolleys and stretchers for all care environments

Manufactured in the UK to the same high standards as the Freeway Shower Chair Range, the comprehensive range of Freeway Shower Stretchers and Trolleys allow carers to shower, change and transfer clients in a safe and dignified manner.

Multiple fixed height and adjustable height options are available throughout the range to ensure carers are working in the safest possible position in order to reduce the risk of injury to themselves and the client.

The Freeway Shower Trolley and Stretcher range has been long established as a market leader due to their durability and versatility, making them suitable for use in a wide range of care environments including hospitals, care homes and individual residences.

All standard models are weight tested to 200kg (30 stone) and come with a standard 12 month warranty.



Shower Stretcher



Freeway Wall Mounted Shower Stretcher

The Freeway Wall Mounted Shower Stretcher is our economy shower stretcher with vinyl top and drainage holes. Photo shows optional padded cot side.

Available in six sizes from 1400mm to 1900mm and can be adjusted on installation up to a working height of 900mm.

Freeway Wall Mounted Shower Stretcher with Padded Roll Edge

This Freeway Wall Mounted Shower Stretcher benefits from a padded roll edge for safety and comfort which often means the cotside is not necessary. Available in six sizes from 1400mm to 1900mm and can be adjusted on installation up to a working height of 900mm. Option to add some sag to the vinyl or fit a skirt.





Freeway Wall Mounted Shower Stretcher with Adjustable Backrest

This Freeway Wall Mounted Shower Stretcher is available with adjustable backrest and can be left or right handed. Available in six sizes.



Wall Mounted Padded Shower Stretcher & Changing Table

This Freeway model serves two purposes in one as it can be used as both a shower stretcher and a changing table. It benefits from a comfortable 50mm vinyl upholstered covering and is available in three sizes: 1500mm, 1700mm and 1900mm.

rs & Trolleys



Freeway Fixed Height Shower Trolley

The Freeway Fixed Height Flat Bed Shower Trolley is our economy model which benefits from top quality sealed bearing castors and the option of cot sides on long sides or all round. Available in six sizes; from 1400mm to 1900mm.

The **Freeway Fixed Height Padded Shower Trolley** has exactly the same specification as above but with padded top included. Available in three sizes; 1500, 1700 and 1900mm.

The **Freeway Fixed Height Shower Trolley with Liner** has cotsides all round and liner included (as shown in photo above). Available in three sizes; 1500mm, 1700mm and 1900mm.

The **Freeway Fixed Height Shower Trolley with Adjustable Backrest** is available with adjustable backrest and can be left or right handed.

FreewayAdjustable Height Shower Trolley

The Freeway Adjustable Height Shower Trolley benefits from optional cot sides and is height adjustable by hydraulic foot pump. Six sizes available from 1400mm to 1900mm.

Freeway Adjustable Height Shower Trolley with Adjustable Backrest

This Freeway Adjustable Height Shower Trolley is a variation of the above with the added inclusion of an adjustable backrest. Six sizes available.





Freeway Adjustable Height Trolley with ABS Top & Liner

This Freeway Adjustable Height Shower Trolley benefits from an ABS top and padded liner with an optional pillow.

It is height adjustable by hydraulic pump to a working height of between 600mm to 925mm. Three sizes available; 1500mm, 1700mm, and 1900mm.


Freeway Adjustable Height Padded Changing Trolley

This Freeway Adjustable Height Shower Trolley benefits from a comfortable one piece padded top.

Photo features optional cot sides all round. Three sizes available; 1500mm, 1700mm and 1900mm.



Shower Stretcher Technical



Dimensions based on a standard 1900mm length Shower Stretcher

- 1 = length = 1900mm
- 2 = length with adjustable headrest = 1975mm
- 3 = width = 670mm
- 4 = width when fitted with standard brackets = 810mm
- 5 = width when fitted with short brackets =760 mm
- 6 = width when fitted with cot side = 835/785mm
- 7 = distance from wall to inside of legs=795mm/745mm
- 8 = distance between bracket hinges = 1100mm
- 9 = height = 900mm

10=distance from wall when folded up=120mm

NB: Padded rolled edge increases width/length by 50mm

Shower Trolley Technical



- $8 = \max$ working height = 925mm
- $9 = \max$ height with cot sides fitted = 1125mm
- 10 =length of trolley base = 1500mm

Fixed Height Shower Trolley

- 2 =length with adjustable headrest = 1975mm
- 3 =length with cot ends fitted both ends = 1950mm

1 and 3

- 5 = width with cot sides both sides = 720mm
- 7 = height with cot sides fitted = 1100mm

Wedged pillow dimensions - width = 550mm, length = 330mm, height = 110mm Waste hose dimensions (for trolley with liner) - length = 1 metres, diameter = 75mm



Contact Us Units 5 and 6 Lower Vicarage Road Southampton Hampshire SO19 7RJ T: 02380 420678 E: info@easycaresystems.co.uk



www.easycaresystems.co.uk