



harvesthealthcare
a prismhealthcare company



Product Code: HWB004

Description:

- The Bed Side Platform Bumper has been designed to easily fit onto the outer edge of the bed frame.
- To reduce the risk of injury for patients on the frame when getting in and out of bed.
- Suitable for use on all our Woburn bed range and similar types of profiling beds produced by other manufacturers.
- Created using soft to the touch material, protecting the user and increasing patient safety.

Cleaning and Care:

The cover is removable and machine washable, for better infection control. Pre wash cold 10 minutes
Main wash 80°C 10 minutes
Followed by cold rinses and extraction.
Do not use abrasive cleaners, phenol disinfectants, solvents, or alcohol-based cleansers, e.g. Dettol, Phenicol, Hibiscrub, Clearsol, Stericol, Hycoline, as these will damage the cover materials.
Tumble dry at low heat.
Do not iron.

Ensure that the cover is thoroughly dried before inserting the foam and placing the bumper into storage.

Latex free: Yes

User Instruction: Yes

Flammability: BS 7175: 1989 (2013) ignition sources 0,1

www.easycare systems.co.uk

Dealer Contact Details:
Easy Care Systems Ltd

Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk