


## 3-LAYER VISCO TOP CUSHION



<b>Product Code:</b>	HFC129.01
<b>Description:</b>	<ul style="list-style-type: none"> <li>• A high risk viscoelastic and foam bonded cushion with improved immersion of visco top and a combination of a second layer of viscoelastic foam and a firm foam base add support.</li> <li>• The foam combination enhances pressure distribution and combines stability and comfort for the patient at high risk.</li> <li>• Two way stretch cover with antislip base.</li> <li>• Cover is latex free and clearly printed, including cleaning instructions.</li> </ul>
<b>Dimensions:</b>	430 x 430 x 75mm
<b>Max. User Weight:</b>	158kg/25 stone
<b>Pressure Relief Risk:</b>	High Risk
<b>Cover: (Available Separately)</b>	Four-way stretch, vapour permeable, PU coated top in navy blue. Zip fastened on 2 sides with single flap cover to prevent fluid ingress. Base section is antislip PVC to reduce unwanted movement.
<b>Foam:</b> Density and Hardness	Top: 58-62kg/m <sup>3</sup> 70-100N Middle: 58-62kg/m <sup>3</sup> 100-130N Base: 34-37kg/m <sup>3</sup> , 90-120N
<b>Fire Retardancy:</b>	BS 5852-1: 1988 crib 5



Visco Top Foam: 58-62kg/m<sup>3</sup> 70-100N

Visco Foam: 58-62kg/m<sup>3</sup> 100-130N

Firm Foam Base: 34-37kg/m<sup>3</sup>, 90-120N

### **Cleaning Guidance:**

Wipe down cover with hot water and mild detergent. Or machine wash and dry. Should the foam interior become damp remove cover and allow to dry naturally. Use no bleach. Dry completely before storage. For a comprehensive list of tested cleaning products see our website.

### **Warranty**

In terms of Harvest Healthcare Guarantee and Warranties the 3-Layer Visco Top Cushion is guaranteed for three years from date of purchase with the following exclusions:

#### **Exclusions:**

1. Removable cushion covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded



#### **Dealer Contact Details:**

Easy Care Systems Ltd

Units 5 and 6 Lower Vicarage Road  
Southampton  
Hampshire  
SO19 7RJ  
T: 02380 420678  
E: [info@easycare systems.co.uk](mailto:info@easycare systems.co.uk)