

## KENSINGTON ACTIVE OVERLAY MATTRESS



**Product Code:**

HAMKENS001

**Description:**

- The Kensington is Harvest Healthcare's Active Overlay Mattress
- 18 cells made of TPU material to reduce noise when the cells inflate and deflate
- Pillow section for patient stability
- All cells are interchangeable for ease of repair and cleaning
- Light mattress enabling the carer to move the mattress easily
- 4-way stretch welded PU cover which is vapour permeable, water resistant and has a Crib 5 fire rating
- Compatible with all of Harvest Healthcare's pumps allowing for a range of maximum patient weights

**Mattress Dimensions:**

2000 x 900 x 110mm

**Weight:**

6kg

**Max. Patient Weight:**

Harvest 1 Pump: 25 Stone (158 kg)  
Harvest 2 Pump: 29 Stone (184kg)

**Pressure Relief Risk:**

High

**Pump Technical Specification:**

Please see separate pump specification sheet

**Cleaning Guidance:**

Cover should be wiped down with hot water and mild detergent. If heavily soiled, the cover may be machine washed up to 65°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

**Warranty**

In terms of Harvest Healthcare Guarantee and Warranties the Harvest Kensington mattress is guaranteed for three years from date of purchase.

**Exclusions:**

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded.



Special Products and Service

**Contact Us**

Units 5 and 6 Lower Vicarage Road  
Southampton  
Hampshire  
SO19 7RJ

T: 02380 420678

E: [info@easycare systems.co.uk](mailto:info@easycare systems.co.uk)

[www.easycare systems.co.uk](http://www.easycare systems.co.uk)