



Product Code: HFB002

Description: Safe Roll Ultra Low- Floor Bumper

- The Safe Roll Ultra Low bumper has been designed to fit onto our Woburn Ultra-Low bed.
 - The product fills the need to provide a safe environment for the user who may not need side rails fitting.
 - The Safe Roll Ultra Low bumper helps to prevent accidental contact with the floor, the material we use is soft to the touch giving the user protection in case they accidentally roll off the bed.
 - Ideal for nursing homes and care in a domestic environment.
 - It is easy to fit. Fitting instruction video is available on the website.
- The cover is removable and machine washable, for better infection control.

Cleaning and Care:

Pre wash cold 10 minutes
Main wash 80°C 10 minutes
Followed by cold rinses and extraction.
Do not use abrasive cleaners, phenol disinfectants, solvents, or alcohol-based cleansers, e.g. Dettol, Phenicol, Hibiscrub, Clearsol, Stericol, Hycoline, as these will damage the cover materials.
Tumble dry at low heat.
Do not iron.
Ensure that the cover is thoroughly dried before inserting the foam and placing the bumper into storage.

Latex free: Yes

User Instruction: Yes

Flammability: BS 7175: 1989 (2013) ignition sources 0,1

Dealer Contact Details:
Easy Care Systems Ltd

Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk