



harvesthealthcare
a prismhealthcare company



Product Code:

HCSBUMP112	Standard 2 Bar Bumper with Net/Mesh Inserts for Woburn, Community & 1200 Bed (L: 2000 mm)
HCSBUMP212	Extended 2 Bar Bumper with net/mesh inserts for Woburn, Community & 1200 Bed (L: 2200 mm)
HCSBUMP130	Standard 3 Bar Bumper with Net/Mesh Inserts for Woburn Ultimate Bed (L: 2000 mm)
HCSBUMP112 + HCSBUMP112.EX20	Extended Standard 2 Bar Bumper with Net/Mesh Inserts for Woburn, Community & 1200 Bed (L: 2200 mm)
HCSBUMP130 + HCSBUMP130.EX20	Extended Standard 3 Bar Bumper with Net/Mesh Inserts for Woburn Ultimate Bed (L: 2200 mm)

Description:

- The Standard Bumper with net/mesh inserts is designed to fit over side rails, providing a safe protection for users in bed.
- The net/mesh inserts are entwined with the bumpers, to improve visibility between each side rail, reducing the feeling of entrapment.
- The bumper fits all Woburn beds.
- All the bumpers come with covers that are removable and machine washable, for better infection control.

Cleaning and Care:

The cover is removable and machine washable, for better infection control.
 Pre wash cold 10 minutes
 Main wash 80°C 10 minutes
 Followed by cold rinses and extraction.
 Do not use abrasive cleaners, phenol disinfectants, solvents, or alcohol-based cleansers, e.g. Dettol, Phenicol, Hibiscrub, Clearsol, Stericol, Hycoline, as these will damage the cover materials.
 Tumble dry at low heat.
 Do not iron.
 Ensure that the cover is thoroughly dried before inserting the foam and placing the bumper into storage.

Latex free: Yes

User Instruction: Yes

Flammability: BS 7175: 1989 (2013) ignition sources 0,1



2 Bar Anti Standard Bumper



3 Bar Standard Bumper



**Extended Standard Bumper with
Net/Mesh Inserts**

Dealer Contact Details:

Easy Care Systems Ltd

Units 5 and 6 Lower Vicarage Road

Southampton

Hampshire

SO19 7RJ

T: 02380 420678

E: info@easycare systems.co.uk

easy  **care**
systems

www.easycare systems.co.uk