



Product Code:

HFC129
HFC123

Description:

- A 50mm layer of viscoelastic foam tops the standard foam base, transforming the cushion to a high risk product.
- The top layer of foam provides extra support and comfort to high risk patients.
- It has a navy blue, four-way stretch cover with an anti-slip base.

Dimensions:

HFC 129: 430x 430 x 76mm
HFC 123: 510 x 510 x 76mm

Weight of Cushion:

HFC 129: 0.8 kg
HFC 123: 1.1 kg

Max. User Weight:

HFC129: 127 kg
HFC123: 139kg

Pressure Relief Risk:

High Risk

**Cover:
(Available Separately)**

Vapour permeable, antibacterial PU coated 4 way stretch cover in navy blue with zip fastening. Silver based antibacterial properties in the PU coating are wash resistant and environmentally friendly.

**Foam:
Density and Hardness**

Base 32-34kg/m³ 155-195N Top 48-52kg/m³, 70-100N

Fire Retardancy:

BS 5852-1: 1988 crib 5

Cleaning Guidance:

Wipe down cover with hot water and mild detergent. Or machine wash and dry. Should the foam interior become damp remove cover and allow to dry naturally. Use no bleach. Dry completely before storage. For a comprehensive list of tested cleaning products see our website.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Visco Top Cushion is guaranteed for three years from date of purchase with the following exclusions:

Exclusions:

1. Removable cushion covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

**Dealer Contact Details:**

Easy Care Systems Ltd

Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk

www.easycare systems.co.uk