

www.easycaresystems.co.uk

## WHEELCHAIR MARINA-ANFI


- **Great buoyancy and stability:** it has two floating armrests and of wide wheels that allow the displacement on any place (beaches, pools with ramp ...)
- \*Comfortable and easy to use: facilitates the transfer of people with reduced mobility for bathing or resting in a very comfortable.
- **Design and comfort:** made of high strength materials, easy to clean, disassembled and folded quickly for transport or storage.
- **Security:** must be used with the assistance of another person, respecting the safety precautions and current regulations.



## MAIN PARTS OF THE WHEELCHAIR



- The backrest and foot rest are mesh polyester, which favors:
  - the outflow of water and sand.
  - greater ventilation and comfort of the person sitting.



- 1. Accompanying handle
- 2. Backrest
- 3. Side floater
- 4. User handle
- 5. Seat
- 6. Foot rest
- 7. Traction bar
- 8. Carry structure
- 9. Wheels

## **TECHNICAL CHARACTERISTICS**

PRODUCT	REFERENCE	SIZES IN MM	WEIGHT OF CHAIR	MAXIMUM WEIGHT OF THE USER	MATERIALS	BACKREST AND SEAT	WHEELS	
MARINA ANFI	102605012	1230×1200×2600	31,7 kg	130 kg	Aluminum Stainless steel	Polyester fiber coated with PVC	Diameter 380 mm.	Inflation pressure 2 kg/cm²
		640×560×1220			Polyethylene Polyurethane			

Conforms to standards: UNE-EN 12182 UNE-EN 12183 Directive 93-42-CEE 2 years warranty





