

Patient Handling Belt


The Patient Handling Belt offers carers a means of providing standing assistance to patients. It can be used independently or in combination with a compatible manual standing/transfer device, e.g. the Oxford Switch. It may also be used to help with walking rehabilitation/gait training activity following careful assessment.



- ✓ Side handles for standing assistance and support
- ✓ Adjustable belt enables use with the Oxford Switch
- ✓ Non-slip inner layer ensures secure fit
- ✓ Elasticated top band for increased comfort
- ✓ Colour coded binding indicates belt size
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed
(From/to seated position)



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SPECIAL	180 kg
S	SL3151	180 kg
M	SL3152	180 kg
L	SL3153	180 kg
XL	SPECIAL	180 kg

