



Repose Cover Solutions

Repose Cover Solutions

Enhanced patient comfort and
added product protection

repose[®]

Repose Covers Solutions

Enhanced patient comfort and added product protection

The Repose range adopts an effective method of pressure redistribution and loading by reducing the pressure exerted at the interface between the patient and the supporting surface.

Repose Covers Solutions add extra protection to Repose products helping to keep Repose in pristine condition.

Repose Companion and Trolley Companion Covers are supplied with a slide sheet base material combines polyurethane and nylon for durability and is welded to the multi stretch cover. Improving patients sense of comfort and security throughout their journey.

Repose Cover Material

Conformable, polyurethane cover material which is durable, latex free, easy to clean and maintain.

Repose covers come with a One-year Manufacturing warranty.

Adjustable Retention Straps

Adjustable securing straps, to ensure Repose products are fitted securing in place aiding stability and security.



Repose Mattress Overlay Cover

Repose Cushion Cover

Repose Care-Sit Cover



Repose Contur Cover

Repose Companion Cover



Safe & Secure Fit



Lateral Transfer

The slide sheet base material combines polyurethane and nylon for durability and is welded to the multi-stretch cover. Improving patient's sense of comfort and security throughout their journey.



Lateral Patient Transfers



Cover Print

Washing instructions and identification boxes can be found on Repose Covers which provides easy visual guidance.

Cover Zips

Conventional cover zip with integral cover flaps or waterproof zip for easy insertion of Repose, easy access to the inflation valve and to protect against fluid ingress.

Cover zips protect against fluid ingress



Multi-patient Use

Repose covers are easy to clean in line with local policy and reusable for multi-patient use, increasing the lifetime of the product.



Multi-Patient Use Product

Easy to clean and reuse





Special Products and Service

Smart valve technology provides reassurance that your Repose product is never over inflated when following the instructions for use.

For more information:

Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk

repose[®]