



Repose[®] Flex

Designed to reduce
the risk of DVT¹

repose[®]

Repose® Flex



repose[®]

Utilising immersion and envelopment, Repose Flex is designed to place the knee in slight flexion (5° to 10°) and allow complete heel offloading while effectively re-distributing the weight of the patient's leg across its surface².

By placing the knee in slight flexion, the risk of Popliteal Vein Compression at the rear of the knee is reduced¹.

Popliteal Vein Compression is a significant source of Deep Vein Thrombosis for patients undergoing surgery or prolonged bed rest³.



Suitable for:

- Patients identified as being at risk of developing a DVT
- Bedbound patients, or patients undergoing surgery
- Protecting patients at all levels of risk of developing all stages of pressure ulcer
- Treatment of pressure ulcers when used as part of a package of care
- Deflatable for convenient storage and transport
- Repose smart valve technology for optimum inflation levels
- Three years warranty as standard



Repose protects patients at all levels of risk



Repose is indicated for the treatment of all categories of pressure ulcers



Multi-Patient Use Product



Repose is supplied with a 3 year warranty as standard

Description	Size (L/W/H)
 Repose Flex & Pump - Extra Small	L 420mm x W 700mm x H 90mm
 Repose Flex & Pump - Small	L 460mm x W 700mm x H 90mm
 Repose Flex & Pump - Standard	L 500mm x W 700mm x H 90mm

For more information:

Units 5 and 6 Lower Vicarage Road Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk

easyCare
systems

Special Products and Service