

Special Products and Service

Ø **Repose Lite Cushion** ŕ

Repose[®] Lite Cushion

Pre-inflated Reactive Air Cushion



Repose[®] Lite Cushion

Repose[®] Lite Cushion is a unique reactive air, single cell, pre-inflated pressure relieving and redistribution cushion with securing strap.

Indicated for the treatment of all categories of pressure ulcers as part of a pressure ulcer prevention strategy, and for patients identified as being very high risk.

The cushion is compatible with all types of static chairs where it can be secured in place with the integrated fixing strap. It can also be used with general purpose wheelchairs with fitted armrests.

The Repose Lite Cushion can be cleaned and reused for multiple patients.

- Pre-inflated
- High levels of immersion and envelopment for effective pressure redistribution and comfort
- Can be used for treatment of pressure ulcers when used as part of a package of care
- Suitable for static chairs and general purpose wheel chairs
- Cost Effective
- Easy to clean
- TPU Strap
- The Repose Blue Cushion Cover available offering an additional level of protection







Repose is indicated for the treatment of all categories of pressure ulcers



One-Year Warranty



Weight Limit 139kg

Multi-Patient Use Product



A Pre-inflated Solution

- Repose Lite Cushion is pre-inflated, ready for use, with a simple one step installation.
- Repose Lite Cushion offers an instant solution for pressure ulcer prevention as part of a package of care.



The Repose Lite Cushion can be used in conjunction with a Repose Blue Cushion Cover offering an additional level of protection.

The Repose Lite Cushion and Repose Blue Cushion Cover are also available packaged together as a complete set.



Special Products and Service

Repose Lite Cushion TPU Strap

The integrated TPU strap secures the product firmly in place aiding stability and security.



Higher Outer Air Cells

Higher outer air cells provide greater stability, additional pelvic support, increasing patient comfort, confidence and positioning.

Low Profile Cells

The low-profile cells prevent patients sitting too high on a chair aiding patient safety and increasing comfort.



easybare systems

Special Products and Service

For more information:

Units 5 and 6 Lower Vicarage Road Southampton Hampshire SO19 7RJ T: 02380 420678 E: info@easycaresystems.co.uk

