

# Welcome®-Giant Half-Step





## Until now, the traditional half-step simply offered an ambulant user shallower step access.

The proven and now ALL-NEW Welcome®-Giant Half-Step reaches a new level!

This totally unique and intelligent range takes into account that people with various walking challenges need more workable solutions.







## 3 sizes...

## MIDI



45cm (18ins)

## **LARGE**



60cm (24ins)

## EXTRA-LARGE



90cm (36ins)

## To suit all requirements:

- Walking stick? If so, the MIDI is perfect!
- Walking frame user? Take a look at the LARGE!
- Is a rollator being used? The EXTRA-LARGE is the answer!

Take a look at page 6 for the easy product selector!

Quick Tip: The bigger the half-step, the safer and easier it is to use!



## With or Without?

Useful Accessory Options - So Easy to Choose and Use!

- Handrails none, one side, or both sides!
- Product kits for easy ordering

## 3 challenges solved...

- 1. Safe and Manageable **Step Heights**
- 2. Eliminate the Trip Risk of Raised Threshold Obstacles
- 3. Instant Safe Handrails

## **Even More User Benefits**

- Safely pause at the door
- All door types including PVC
- No building work required

## 1. Safe and Manageable **Step Heights**



A single Welcome®-Giant Half-Step will instantly create a safe and manageable step for thresholds up to 30cm (12ins). Add handrails and other accessories as required for the specific user and application.



Unlike any other, the Welcome®-Giant Half-Step can be installed level with the top of the threshold and an optional doorplate accessory used to create a flush transition with no unsafe gaps or 'step-over' obstacle risk. This especially improves access for those with limited leg/knee mobility, drop foot or shuffling gait. The user simply steps up or steps down, rather than risking negotiating up and over the raised threshold.

## 3. Instant Safe Handrails



The multiple benefits of the Welcome\*-Giant Half-Step solution include providing safe and instant handrails wherever required. Perfect for situations when attaching permanent grab rails or handrails to a property is not practical. It also removes the risk of falls from leaning over the step to reach a wall fixed grab rail.

The doorplate accessory is a great solution for curved or rounded steps.

## **LARGE MIDI** Easy product selector 75cm (30ins) 75cm (30ins) 75cm (30ins) 45cm (18ins) 60cm (24ins) 90cm (36ins) **Ambulant** All sizes are suited for ambulant use. Limited Knee/Leg **Mobility** Larger sizes are more suited for limited knee mobility. **Walking Frame** Larger sizes suited to most walking frames. Rollator Extra-Large size is especially suited to most rollator devices. **Limited Space** Where space is limited, choose Large or Midi sizes. **Outward Opening Doors** Extra-Large half-step provides more space for safely pausing to open the door. Capacity The high capacity is also perfectly suited for bariatric use. Important note: 400kg capacity requires all feet to be securely fixed to the ground 400kg 400kg 400kg and requires connected handrails both (63st) (63st) (63st) sides if using height extension kit WHSLXK. Capacity is otherwise 300kg (47st).

Specifications and Accessories HALF-STEPS		MIDI	LARGE	EXTRA- LARGE
ACCESS	ORIES			
	<b>Handrail</b> With integral safety kerbs. Fit to any side or opposite sides. Universal height: 90cm (36ins)	<b>WHS45</b> 45cm (18ins)	<b>WHS60</b> 60cm (24ins)	WHS75 75cm (30ins WHS90 90cm (36ins
	90cm Anti-Trin Too Plate			



#### 90cm Anti-Trip Toe Plate

Provides additional trip protection on the long (90cm) Extra-Large Half-Step sides WHSTP90



### Doorplate

Closes the gap between the Giant Half-Step and the door frame Self-supporting for PVC frames WHSDP75

75cm (30ins)



#### **Height Extension Kit**

Extend any Giant Half-Step up to 20cm (8ins) high. Useful for 'building over' existing steps

**WHSLXK** 



### **Handrail Joiner Kit**

Join handrails of any two or more adjacent half-steps (2 pieces)

**WHSRJK** 

KITS FOR EASY ORDERING		MIDI	LARGE	LARGE
	<b>Giant Half-Step Kit</b> Giant Half-Step Kit with Hand Rail. Height range: 8 - 15cm (3.1 - 6ins)	WHSMR Includes: WHSM WHS45	WHSLR Includes: WHSL WHS60	WHSXLR Includes: WHSXL WHS90
	<b>Double Giant Half-Step Kit</b> Double Giant Half-Step Kit with Hand Rail. Height up to: 20cm (8ins)	WHSM2R Includes: 2 x WHSM 2 x WHS45 WHSRJK WHSLXK	WHSL2R Includes: 2 x WHSL 2 x WHS60 WHSRJK WHSLXK	WHSXL2R Includes: 2 x WHSXL 2 x WHS90 WHSRJK WHSLXK





Special products and Service

## **Contact Us**

Units 5 and 6 Lower Vicarage Road Southampton, Hampshire, SO19 7RJ T: 02380 420678

E: info@easycaresystems.co.uk