

## Argyll™ II with CLP

The Argyll™ II is a high-end, specialised dynamic mattress replacement system with a digital pump and CLP (Constant Low Pressure), targeted towards assistive treatment for patients who have existing and potentially severe pressure damage.


This product was initially designed for very high risk patients in an acute setting. However, the product is well adapted for use in long term care.

The 8 inch (200mm) high cells, cyclically inflate and deflate in an alternating pattern, providing direct pressure relief throughout the dynamic process. Additionally, the mattress incorporates CLP (Constant Low Pressure) as an alternative therapeutic support surface which maintains stability across the surface area at all times.

**If you have any questions we are here to help.  
Call us on 0800 193 6030.**

### PRODUCT FEATURES

- Designed for Very High Risk patients
- Suitable for both acute and long term care settings
- Dynamic inflation/ deflation cycle provides direct pressure relief
- Incorporates Constant Low Pressure (CLP) as an alternative therapeutic supports surface
- Cell on cell construction base cells stay permanently inflated in case of power cut
- Dual mode: alternating/CLP



## Features and Benefits

- Multi-stretch, vapour permeable, antimicrobial, cover material aids infection control
- Fully High Frequency welded product prevents fluid ingress
- Large/Deeper Cell design aids patient comfort, whilst providing excellent pressure relief
- CLP (Constant Low Pressure) allows patient immersion on a non-moving surface.
- CPR Function aids rapid deflation or medical attention
- Transport Mode. 48 hour transport function. Enables ward/patient procedures.
- UK Manufactured product (cover) provides total confidence to user via exceptional quality and workmanship
- Clinically proven, Giving peace of mind to nursing staff and patients alike
- Static Mode aids patient transfers and reverts back after 30 minutes

## TECHNICAL INFORMATION

Product Code	ARG-142-2 For dual mode with CLP use code ARG-142-2-CLP
Risk Category	Very High Risk
Fatigue Class	Very Severe
Maximum Patient Weight	35 stone 222kg
Alternation cycle	1 in 2
Alternation time	11.5 minutes (1 in 2)
Alternation Cycle	1 in 2
Alternation Time	11.5 mins
Mattress Dimensions	200cm x 88cm x 20cm (suitable for most hospital/homecare beds)
Mattress Weight	Less than 7kg
Mattress Cover Coating	200-230g
Warranty	2 years – extended warranty available

WIN  CARE

easy  care  
systems

Special Products and Service

Contact Us  
Units 5 and 6 Lower Vicarage Road  
Southampton  
Hampshire  
SO19 7RJ  
T: 02380 420678  
E: [info@easycare systems.co.uk](mailto:info@easycare systems.co.uk)