



Product Codes:

HSMCOM
HSMCOM.12

Description:

- The Community Mattress is a cost effective comfort static mattress for use on profiling beds.
- The high resilient, fire retardant foam is covered with a zipped four-way stretch, water resistant PU cover and is shaped to allow profiling.
- The community mattress provides excellent pressure distribution and can be turned to extend service life.

Mattress Dimensions:

2000 x 900 x 125mm
2000 x 1200 x 125mm

Mattress Weight:

8kg
10kg

Max. User Weight: Risk

25 Stone (159kg)

Level:

Medium

Foam Density Range:

31-34 (kg/m³)

Foam Hardness Range:

100-130N

Fire Retardancy:

BS 7177 Ignition Source 5

Cover:

(Also available separately)

Four-way stretch PU coated polyester cover, with zip fastening. The coating is water resistant and antibacterial as well as being vapour permeable, washable and flame retardant (crib 5).

Cleaning Guidance

The mattress foam may be hand washed with mild detergent. Do not immerse in water. Allow to dry naturally. Cover should be wiped down with hot water and mild detergent or machine washed up to 95°C and tumble dried. For a comprehensive list of tested cleaning products see our website.

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Community mattress is guaranteed for one year from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded



Special Products and Service

Contact Us

Units 5 and 6 Lower Vicarage
Road

Southampton
Hampshire

SO19 7RJ

T: 02380 420678

E: info@easycare systems.co.uk

www.easycare systems.co.uk