



WOBURN COMMUNITY LOW PROFILING BED



Product Code:	HLBC791L HLBC791LAR (with Auto-Regression)
Description:	<ul style="list-style-type: none"> • Built to the same specification as our popular Woburn Community Bed • Offers the benefit of a lower entry and exit height • Ideal for many environments, fully compliant with legislation. • Comfort adjustments of height, back rest or knee break done simply using the lockable handset.
Safe Working Load:	220kg / 34 Stone
Max. Patient Weight:	185kg / 29 Stone
Overall Dimension:	2115 x 1015mm
Mattress Base:	2000 x 900mm
Max. Mattress Depth:	205mm
Adjustable Height:	220mm-640mm
Angle of Back Rest:	Up to 70 degrees
Angle of Knee Break:	Up to 30 degrees
Bed Weight:	81kg + 11kg with rails
Key Information:	4 Section Profiling Bed Electronically Powered Back Rest & Thigh Rest Manual Leg Rest Adjustment Wooden Head and Foot Boards Side Rails - Optional Auto-Regression and Inbuilt Length Extension - Optional Trendelenburg Function as standard Lockable Handset Low Voltage – 29v; 9v Battery Backup Comprehensive Manual Warranty: 5yrs

Supplier Contact Details:

t: +44(0)1709 377172

e: enquiries@harvesthealthcare.co.uk

w: www.harvesthealthcare.co.uk

Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60



Special Products and Service

Contact Us

Units 5 and 6 Lower Vicarage Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk

www.easycare systems.co.uk