

1)





Repose Contur Riser recliner pressure redistribution overlay



## **Repose Contur**

## Repose Contur: pressure area care with enhanced comfort

It has long been recognised there is a need for a clinical solution to the problem of pressure damage through the extended use of riser recliner chairs by at-risk patients.

Even when a pressure redistribution cushion is built into a chair, this does not address the potential damage that may be caused to other vulnerable parts of the patient's skin such as the scapulae, occiput or heels.

Repose Contur is a pressure redistribution cushion overlay for riser recliner chairs, utilising the tried and trusted Repose inflatable air cell and 'smart valve' technologies.

Offering pressure area care to the back, seat and leg sections, for those who require a specialist support surface. Repose Contur secures to all type of chair sizes easily with retention straps. Providing patient security and comfort.



Repose products protect patients at all levels of risk

VERY HIGH RISK 



Repose is indicated for the treatment of all categories of pressure ulcers



The Repose inflation pump with Smart Valve Technology, ensures Repose products are

inflated to the correct pressure offering high

levels immersion and envelopment.









Special Products and Service

## **Contur Retention Straps**

Adjustable and secure strap system allows for a secure fit on a range of chair sizes



Three Section Design

Conforms to the three sections of a riser recliner chair. Enhancing patient comfort, offering top to toe pressure redistribution.

Repose Contur Cover

Conformable polyurethane cover which is durable, easy to clean and maintain.



Repose Contur Weight Limit 139kg

## easybare systems

Special Products and Service

For more information:

Units 5 and 6 Lower Vicarage Road Southampton Hampshire SO19 7RJ T: 02380 420678 E: info@easycaresystems.co.uk

