



Repose[®] Cushion
Comfortable pressure
redistribution

repose[®]

Repose® Cushion

Comfortable pressure redistribution

Indicated for the treatment of all categories of pressure ulcers as part of a pressure ulcer prevention strategy, and, for patients identified as being very high risk.

The cushion is compatible with all types of chairs where securing in place with the strap is possible and can also be used with general purpose wheelchairs with fitted armrests.

The Repose Cushion can be cleaned and reused for multiple patients.

- Portable solution for static chairs
- Slim profile increases patient comfort
- Higher outer air cells provides greater stability
- Three year manufacturing warranty
- Repose Pump with Smart Valve included
- Easy to clean for multi-patient use



Repose products protect patients at all levels of risk



Repose is indicated for the treatment of all categories of pressure ulcers



Multi-Patient Use Product



Weight Limit 139kg



Repose is supplied with a 3 year warranty as standard



Repose Cushion Retention Strap



Securing the product firmly in place aiding stability and security.



Smart Valve



Repose Pump & Smart Valve

The Repose pressure redistribution cushion comes packed inside a unique and easy to use pump, which enables the inflatable air cushion to be inflated ready for use within seconds and automatically ensures that the product is inflated to the correct pressure.

Higher Outer Air Cells

Higher outer air cells provides greater stability, increasing patient comfort, confidence and positioning.

Low Profile Cells

Patient doesn't sit too high on the chair aiding patient safety and increasing comfort.

Accessories

The Repose Cushion has the following accessories; Repose Refresh Cushion Cover, Repose Strap Extender and the Repose Blue Cushion Cover.

All accessories have a one year warranty from date of despatch.



Blue Cushion Cover



Safe & Secure Fit





Special Products and Service

For more information:

Units 5 and 6 Lower Vicarage
Road
Southampton
Hampshire
SO19 7RJ
T: 02380 420678
E: info@easycare systems.co.uk

repose[®]